



## FEATURES AND BENEFITS

- Preloaded fully disposable delivery system with fully contained sutures
- Two piece 'Stable Drive' handle design
- PEEK and Titanium anchors available in 5.5mm and 6.5mm diameters
- Integrated stainless steel punch for insertion without pilot hole preparation
  - Stainless steel punch removed upon implant release
- Double lead tapered threads for faster delivery
- Three #2 High Strength UHMWPE/Polyester sutures

## ORDER INFORMATION

ITEM #	DESCRIPTION
HNM-78-5515	5.5mmX15mm
HNM-78-6515	6.5mmX15mm
HNM-91-5515-Ti	5.5mmX15mm, Ti
HNM-91-6515-Ti	6.5mmX15mm, Ti

