



CANNULA SYSTEM

The HNM Cannulas are available in varying diameters and lengths for arthroscopic shoulder and hip surgery. The Cannula's triple-spaced flex seal technology is designed to hold water in the joint, and allows smooth instrument transition.



FEATURES AND BENEFITS

- Thread depth ensures secure tissue retention in the joint
- Color coded obturators for easy differentiation between sizes and visualization inside the joint
 - 6.5mm - Green
 - 7.5mm - Blue
 - 8.5mm - Yellow
- Large ergonomic obturator handles allow simple insertion
- Transparent Cannula body allows easy instrument and suture visualization
- Double lead fully-threaded body design is gentle on tissue
- Sterile Peel-Packed with 6 per box
- Disposable
- 6.5mm, 7.5mm and 8.5mm inner diameter options
- 70mm and 90mm length options



HNM-17-100-6570
Disposable Cannula
6.5mm I.D. x 70mm
Box Quantity: 6



HNM-17-100-6590
Disposable Cannula
6.5mm I.D. x 90mm
Box Quantity: 6



HNM-17-100-7570
Disposable Cannula
7.5mm I.D. x 70mm
Box Quantity: 6



HNM-17-100-7590
Disposable Cannula
7.5mm I.D. x 90mm
Box Quantity: 6



HNM-17-100-8570
Disposable Cannula
8.5mm I.D. x 70mm
Box Quantity: 6



HNM-17-100-8590
Disposable Cannula
8.5mm I.D. x 90mm
Box Quantity: 6