

CARINA MEMORY STAPLE

THE MEMORY
STAPLE IS AN IDEAL
FORM OF FIXATION
FOR ARTHRODESIS,
OSTEOTOMY
AND FRACTURE
FIXATION.




HNM Total Recon is committed in providing surgeons and health care professionals with advanced surgical solutions for effective and time-saving treatment of their patients. Our most important task is meeting surgeons' needs we accomplish this by creating a partnership and driving innovative product development.



Manufactured from nitinol memory alloy, the staple offers internal fixation with residual compression. The wide variety of sizes and styles make it ideal for most applications in the foot, ankle, hand and wrist.

FEATURES & BENEFITS

- Body temperature activated for constant compression.
- S-bend bridge allows for more uniform compression across the entire staple.
- 17 standard sizes.
- 6 Step-staple sizes with offsets of 6mm, 9mm and 12mm.
- Made from Nitinol memory alloy, the Bio-pro memory staple is fully active at 98.6° F (37°C), patient Body temperature
- Simple instrumentation.



SIMPLE INSTRUMENTATION

The Memory Staple's unique activation method is complemented by simple, effective instrumentation that provide the surgeon with several intra-operative options for quick, reproducible results.



Description	Item #
X-Small Staples- 1.2mm Wire Size	
7mmx5mm	HTR-01-769037
9mmx7mm	HTR-01-769138
Small Staples - 1.5mm Wire Size	
7mmx5mm	HTR-01-769328
9mmx7mm	HTR-01-769429
11mm x 8mm	HTR-01-769530
11mm x 10mm	HTR-01-769631
13mm x 10mm	HTR-01-769732
15mm x 12mm	HTR-01-769833
11 mm x 15mm x 13mm	HTR-01-769934
11 mm x 17mm x 15mm	HTR-01-769135
11 mm x 19mm x 17mm	HTR-01-769236
Medium Staples- 2.0mm Wire Size	
15mm x 12mm	HTR-01-769425
18mm x 18mm x 15mm	HTR-01-769526
20mmx20mm	HTR-01-769627
Large Staples- 2.0mm x 3.0mm Wire Size	
20mmx20mm	HTR-01-769822
25mmx22mm	HTR-01-769923
30mm x 30mm	HTR-01-799024
Calcaneal Step Staple- 2.0mm x 3.00mm Wire Size	
15mm x 20mm - 6mm Step	HTR-01-779215
15mm x 20mm - 9mm Step	HTR-01-779316
15mm x 20mm - 12mm Step	HTR-01-779417
20mm x 20mm - 6mm Step	HTR-01-779518
20mm x20mm - 9mm Step	HTR-01-779619
20mm x 20mm - 12mm Step	HTR-01-779720
Instruments	
Drill Bit - 1.4mm, XS	HTR-01-769259
Drill Bit - 1.8mm, SM	HTR-01-769360
Drill Bit, 2.4mm, MD	HTR-01-769761
Drill Bit, 3mm - LG	HTR-01-769162
Anchor Pin	HTR-01-769947
Staple Bulldog Hemo - 4 3/4, SM	HTR-01-769351
Staple Clamp - LG	HTR-01-769452
Staple Punch - LG	HTR-01-769553
Staple Punch - SM	HTR-01-769654
Adjustable Drill Guide	HTR-01-769981
Tray Basic Instrument Set	HTR-01-769008
Tray Silicone Brackets Basic Instrument Set	HTR-01-009100

