



SERPENS
Modular CMC
Implant System

PRODUCT INFORMATION





Contents

Introduction	3
Features	3
Instrumentation	4
Ordering Information	6



HNM Total Recon is committed in providing surgeons and health care professionals with advanced surgical solutions for effective and time-saving treatment of their patients. Our most important task is meeting surgeons' needs. We accomplish this by creating key partnerships and driving innovative product development.

📍 20855 Northeast 16th Avenue, Suite C15
Miami, Florida 33179

☎ 1.866.291.8498 ✉ info@hnmtotal.com
📞 1.888.608.9788 🌐 www.hnmtotal.com


SERPENS MODULAR CMC

IMPLANT SYSTEM

The SERPENS Modular CMC Implant provides relief for arthritis in the metacarpal-trapezium joint. Arthritis of the thumb is very common and can be severely disabling. Relief of pain and improved grip strength can be expected with the SERPENS Modular CMC Implant. The design of the medial offset head and varus angle of the implant stem provide excellent stability and reduce the likelihood of dislocation and offers greater range of motion than other implant options. Additionally, the modular design offers the surgeon intra-operative flexibility to individually fit the stem to the metacarpal and the head to the trapezium, while maintaining natural joint tension.

The SERPENS Modular CMC Implant has a varus angled stem which helps prevent dislocation. Also the titanium plasma sprayed stem allows for bony ingrowth. The implant restores range of motion and provides pain relief, with little to no effect on cosmetics of the hand.

FEATURES







INSTRUMENTATION

SERPENS MODULAR CMC BASIC INSTRUMENT SET

HTR-01-000666



CONTENT

BASIC INSTRUMENTS

Item #	Description	Qty	Item #	Description	Qty
HTR-01-520359	Modular Thumb Impactor	1	HTR-01-940750	Allen Key 2.00mm & 5/32	1
HTR-01-720260	Trial Head 12mm	1	HTR-01-940861	Stainless Steel Wrench 1/4	1
HTR-01-720361	Trial Head 13mm	1	HTR-01-940962	Silicone HDL (4 3/4 SQC Cann 1.8mm)	1
HTR-01-720462	Trial Head 14mm	1	HTR-01-941063	Cutting Guide	1
HTR-01-720563	Trial Head 15mm	1	HTR-01-950132	Impactor 7.5 mm	1
HTR-01-720664	Trial Head 12mm+2	1	HTR-01-950233	Impactor 8.5mm	1
HTR-01-720765	Trial Head 13mm +2	1	HTR-01-950334	Impactor 10.0 mm	1
HTR-01-720866	Trial Head 14mm+2	1	HTR-01-950435	Impactor 11.5mm	1
HTR-01-720967	Trial Head 15mm+2	1	HTR-01-970558	Burr Course 12mm	1
HTR-01-761013	Impactor Tip (curved)	1	HTR-01-970659	Burr Course 13mm	1
HTR-01-760114	Trial Stem 7.5mm	1	HTR-01-970760	Burr Course 14mm	1
HTR-01-760215	Trial Stem 8.5 mm	1	HTR-01-970861	Burr Course 15mm	1
HTR-01-760316	Trial Stem 10.0mm	1	HTR-01-980902	Burr Fine 12mm	1
HTR-01-730445	Trial Stem 11.5mm	1	HTR-01-981006	Burr Fine 13mm	1
HTR-01-730546	Burr Adapter - Microaire	1	HTR-01-980107	Burr Fine 14mm	1
HTR-01-730647	Burr Adapter- Stryker	1	HTR-01-980208	Burr Fine 15mm	1
HTR-01-730748	Assembly Block	1	HTR-01-980309	Alignment Guide w/pin	1
HTR-01-730851	Bone Spatula	1	HTR-01-890421	Sizing Guide 12mm	1
HTR-01-730952	Impactor Tip (straight)	1	HTR-01-890529	Sizing Guide 13mm	1
HTR-01-821072	Trial Head 12mm+4	1	HTR-01-910625	Sizing Guide 14mm	1
HTR-01-880117	Trial Head 13mm+4	1	HTR-01-930798	Sizing Guide 15mm	1
HTR-01-880281	Trial Head 14mm +4	1	HTR-01-980820	QC Thumb Trial Impactor	1
HTR-01-730374	Trial Head 15mm+4	1	HTR-01-980932	Head Remover	1
HTR-01-730483	Burr Adapter Synthes	1	HTR-01-001065	Tray for SERPENS Modular CMC	1
HTR-01-730587	Burr Adapter Microchoice	1			
HTR-01-730640	Burr Adapter Mini Microaire	1			



ORDERING INFORMATION

IMPLANT STEM



Description	Item #
7.5mm	HTR-01-751096
8.5mm	HTR-01-750197
10mm	HTR-01-750298
11.5mm	HTR-01-750399

IMPLANT HEAD



Description	Item #
12mm	HTR-01-710499
12mm +2	HTR-01-720538
12mm +4	HTR-01-750600
13mm	HTR-01-700705
13mm +2	HTR-01-720839
13mm +4	HTR-01-750901
14mm	HTR-01-702006
14mm +2	HTR-01-720140
14mm +4	HTR-01-750207
15mm	HTR-01-700307
15mm +2	HTR-01-720441
15mm +4	HTR-01-750508

HNMTOTAL

TOTAL RECON

📍 20855 NE 16 Ave Suite C15, Miami, FL 33179

☎ 1.866.291.8498

📠 1.888.608.9788

✉ info@hnmtotal.com

🌐 www.hnmtotal.com