

ANU VENDAJE





↗ 20855 NE 16 Ave Suite C15, Miami, FL 33179

☎ 1.866.291.8498 ✉ info@hnmtotal.com
📞 1.888.608.9788 🌐 www.hnmmedical.com



ANU LIFE SCIENCES



Amniotic Membrane

Amnion Chorion Membrane

Superiority of placental tissue allograft membranes placental extracellular matrix (ECM) contain:

- Lamin
- Proteoglycans
- Growth factors
- Hyaluronic acid
- Fibronectin
- Collagens 1,2,3,4,5,6,7 all in the bio scaffolds of the ECM, providing an environment for healthy well organized tissue to develop minimizing adhesions, fibrosis and scarring.

The presence of the above components provides an anti-inflammatory environment reducing inflammatory and reducing risk of infection.

It down regulates TGF- β , suppresses pro-inflammatory interleukins, and inhibiting Matrix Metalloproteinases (MMPs)

Anti-Microbial

The presence of integrins and defensins create an environment that potentially can reduce the bacterial colonization of wounds, and surgical sites.

Amniotic Membrane Product description

Biologic allografts derived from dehydrated human amniotic tissue. These products are sterilized post-packaging via e-beam sterilization.



Two thicknesses are available:

- Vendaje amnion membrane approximately 150 microns in thickness
- Vendaje Amnion, Whartons jelly, and chorion approximate-ly 250 microns in thickness

Easy to apply either topically and intra-operatively
Has been shown to minimize scar, adhesions and promote healthy tissue regeneration.

Stored at room temperatures with 5 year shelf life.

Adheres well, no glue, or suture necessary, hydrostatic tension holds in place.

Ideal handling capabilities, provided in dry format, to facilitate placement, then wet, molds and holds position.

Size	Amnion	Amnion Chorion
2x2 cm	200-004-0202	300-004-0202
2x4 cm	200-008-0204	300-008-0204
4x4 cm	200-016-0404	300-016-0404
4x8 cm	200-032-0408	300-024-0406