



# ETHOS BUTTON AL SELF-LATCHING





ADJUSTABLE  
BUTTON  
LOOP

**HNM**   
SPORTS MEDICINE

 20855 Northeast 16<sup>th</sup> Avenue, Suite C15  
Miami, FL 33179

 1.866.291.8498 - 1.888.608.9788

 [www.hnmmedical.com](http://www.hnmmedical.com)

 [Info@hnmmedical.com](mailto:Info@hnmmedical.com)



# ETHOS BUTTON AL


## SELF-LATCHING ADJUSTABLE BUTTON LOOP

### 3-POINT LOCKING DESIGN

Unique 3-point locking design and braid construction for easy placement and to reduce suture creep

### VISIBILITY

Crossover pattern to increase visibility looped inside shiny titanium button



### ONE SIZE FITS ALL

One size fits all your ACL repairs. Self-adjusting loop covering 15 to 60mm ACL loop repair diameters

### OPTIMIZED PACKAGING

Optimized packaging using a single, long, protective card to reduce coiling of the loop and suture

PRODUCT CODE	PRODUCT DESCRIPTION
HNM-OBAL-80-Ti	Ethos Button Al (Adjustable Loop), (2) #5 Suture

📍 20855 NE 16 Ave Suite C15, Miami, FL 33179

☎ 1.866.291.8498

📞 1.888.608.9788

✉ info@hnmmedical.com

🌐 www.hnmmedical.com

**HNM**   
**SPORTS MEDICINE**